


# FUNdamentals of Hematology:

## What, When and Why?

Eric Morissette DVM, Dipl. ACVP (Clinical Pathology)

### A comprehensive complete blood count (CBC) test consists of 2 components

A quantitative CBC and a qualitative blood smear<sup>1</sup>



### Ideally, a blood smear evaluation should be performed as a part of every CBC<sup>1</sup>


At a minimum, blood smears must be performed:

- On every sick patient
- In each instance of abnormal counts or automated cell count flags

Automated cell count flag	Abnormality
Red blood cells (RBCs)	Anemia <sup>2,3</sup>
White blood cells (WBCs)	Cancer, infection, inflammation <sup>2,3</sup>
Platelets (PLTs)	Disease and clumping <sup>3</sup>

### Why aren't blood smears performed more often?<sup>4</sup>

- Lack of experience preparing blood smears
- Time- and labor-intensive process
- Lack of confidence and experience with interpretation
- Assumption that automated count is correct



High-resolution image from VETSCAN IMAGYST™.

A blood smear evaluation **should not be utilized** as a replacement for an automated cell count. If properly maintained, automated analyzers are more precise and accurate than manual cell counts.<sup>5</sup>

# Failure to perform blood smears can lead to errors in clinical decisions

## Blood smears inform clinical decisions and enable veterinarians to<sup>6-9</sup>:

- Confirm automated CBC results
- Assure quality
- Provide additional insights on cell morphology to guide diagnosis and treatment


## Morphological changes that may be identified by a blood smear<sup>1,5,10\*</sup>:

Red blood cells (RBCs)	White blood cells (WBCs)	Platelets (PLTs)
Polychromasia <sup>†</sup>	Left shift (increased neutrophil band cells)	Macroplatelets <sup>†</sup>
Anisocytosis	Toxic changes	
Spherocytes	Reactive lymphocytes	PLT clumping <sup>†</sup>
Heinz bodies	Blast cells	
Fragmented RBCs		
Nucleated RBCs <sup>†</sup>	Mast cells	
RBC parasites		

\*Table includes common examples and is not intended to be an exhaustive list.


<sup>†</sup>Indicates morphological changes currently identified by VETSCAN IMAGYST™ artificial intelligence (AI) blood smear analysis. Other morphology can be assessed via VETSCAN IMAGYST digital cytology image transfer.

## Integrating VETSCAN IMAGYST into a complete, in-hospital hematology solution




**Use any point-of-care hematology analyzer**

The VETSCAN® HM5 is an easy-to-use option that reports a full, 5-part CBC differential with 22 parameters in <4 minutes



**VETSCAN IMAGYST artificial intelligence (AI) technology can review blood smears automatically and quickly**


- Confirm automated cell counts
- Follow up on abnormal automated CBC results
- If abnormalities are observed, expert review via digital image transfer is available\*



**Access expert review by a Zoetis clinical pathologist when needed\***

Digitally submit images for further evaluation beyond AI review, including:

- WBCs—left shifts, toxic changes, malignancy
- RBCs—morphology, inclusions
- PLTs—thrombocytopenia



**Optional complimentary consult**

Obtain free consultations from veterinary specialists with the Zoetis Global Consultation Service whenever further guidance is needed<sup>†</sup>

\*Option to send physical slide to our network of clinical pathologists as needed. Additional costs may apply.

<sup>†</sup>Service available through ZoetisDx platform. Speak to your Zoetis representative to learn more.

**Request a demo today!**



**Learn more**

**References:** 1. Villiers E. Introduction to haematology. In: Villiers E, Ristic J, eds. *BSAVA Manual of Canine and Feline Clinical Pathology*. 3rd ed. British Small Animal Veterinary Association; 2016:27-37. 2. Kahn CM, Line S, Aiello SE. Diagnostic procedures for the private practice laboratory. In: Kahn CM, Line S, Aiello SE, eds. *The Merck Veterinary Manual*. 10th ed. Merck & Co., Inc.; 2010:1487-1492. 3. Barger AM. The complete blood cell count: a powerful diagnostic tool. *Vet Clin North Am Small Anim Pract*. 2003;33(6):1207-1222. doi:10.1016/s0195-5616(03)00100-1. 4. Data on file, VETSCAN IMAGYST Global Acquisition Study, 2021, Zoetis Inc. 5. Harvey JW. Hematology procedures. In: Harvey JW, ed. *Veterinary Hematology: A Diagnostic Guide and Color Atlas*. Elsevier Inc; 2012:11-32. 6. Zabolotzky SM, Walker DB. Peripheral blood smears. In: Cowell R, Valenciano Amy, eds. *Cowell and Tyler's Diagnostic Cytology and Hematology of the Dog and Cat*. 5th ed. Elsevier Inc.; 2020:438-467. 7. Weiss DJ, Tvedten H. The complete blood count, bone marrow examination, and blood banking: general comments and selected techniques. In: Willard MD, Tvedten H, eds. *Small Animal Clinical Diagnosis by Laboratory Methods*. 5th ed. Elsevier Inc.; 2012:12-37. 8. Stirn M, Moritz A, Bauer N. Rate of manual leukocyte differentials in dog, cat and horse blood samples using ADVIA 120 cytograms. *BMC Vet Res*. 2014;10:125. doi:10.1186/1746-6148-10-125. 9. Sharkey L, Heinrich D. In-clinic hematology: the blood film review. *Today's Veterinary Practice*. Published 2015. Accessed January 5, 2022. <https://todaysveterinarypractice.com/in-clinic-hematology-the-bloodfilm-review/>. 10. Weiser G. Laboratory technology for veterinary medicine. In: Thrall MA, Weiser G, Allison RW, Campbell TW, eds. *Veterinary Hematology and Clinical Chemistry*. 2nd ed. John Wiley & Sons, Inc.; 2012:3-33.

All trademarks are the property of Zoetis Services LLC or a related company or a licensor unless otherwise noted.  
© 2022 Zoetis Services LLC. All rights reserved. DIA-00501

**zoetis**